TheraBond[®]3D

antimicrobial barrier systems

United States Distribution Instructions for Use

Description

TheraBond[®]3D Antimicrobial Systems are absorbent, post-operative dressings consisting of silver plated nylon 3-D fabric. TheraBond[®]3D provides an effective barrier to bacterial penetration for up to 14 days. In the presence of exudate the dressing will help maintain a moist wound environment.

Indications

TheraBond[®]3D is indicated for use in light to moderately exuding partial and full thickness wounds including traumatic wounds, surgical wounds, donor sites, 1st and 2nd degree burns, as well as decubitus ulcers, diabetic ulcers and vascular ulcers. TheraBond[®]3D may be used over debrided and partial thickness wounds.

Contraindications

- Do not use on patients with a known sensitivity to silver.
- Do not use on patients during MRI (Magnetic Resonance Imaging) examination.
- Prior to commencing radiation therapy remove TheraBond[®]3D. Apply a new TheraBond[®]3D Antimicrobial Barrier System following treatment.

Precautions

- For external use only.
- TheraBond[®]3D is not compatible with oil base products such as petrolatum.
- TheraBond[®]3D may not be compatible with topical antimicrobials.
- Do not use TheraBond[®]3D with saline solutions or povidone iodine.
- Avoid contact with electrodes or conductive gels during electronic measurements e.g. EEG and ECG.
- TheraBond[®]3D is not intended to provide treatment for infected wounds.
- TheraBond[®]3D may be used on infected wounds which are being managed in accordance with institutional clinical protocols for infection abatement as an adjunct to the standard treatment regimen to provide a barrier to bacterial infection.
- If reddening or sensitization of the skin occurs, discontinue use.
- TheraBond[®]3D should only be used in premature infants (≤ 37 weeks gestation) when clinical benefit outweighs any potential patient risks.

©2015 Alliqua BioMedical, Inc. All rights reserved.

The Alliqua BioMedical logo mark is a trademark of Alliqua BioMedical, Inc.

TheraBond is a registered trademark of Alliqua BioMedical, Inc.

PATENTS PENDING

Manufacturer

Sterile unless package is damaged or open

Sterilization using ethylene oxide

Do not Re-sterilize

Prior to use, see instructions

Latex-free

Do not reuse

Rx Only

Guidelines for Use, Island Systems

Cleanse wound with sterile water.

Choose a System where the pad completely covers the wound.

Remove System from sterile pouch, using aseptic technique.

To apply: remove half of the release paper and place silver pad in contact with the wound.

Remove second half of the release paper and smooth System down.

To remove: lift one corner and slowly remove in a motion parallel to the skin.

Dressing should be changed based on the amount of exudate or the condition of the wound.

Note: The System may adhere on wounds that are dry or have very light exudate. If the System is not easily removed, then moisten or soak to assist removal.

Guidelines for Use, Wraps and Contact Systems

Cleanse wound with sterile water.

Choose a System that completely covers the wound.

Remove System from sterile pouch, using aseptic technique.

To apply: moisten the system with sterile water and place the textured side in contact with the wound. Cover and secure the System per wound care protocol.

To remove: lift one corner and slowly remove in a motion parallel to the skin.

Dressing should be changed based on the amount of exudate or the condition of the wound.

Note: The System may adhere on wounds that are dry or have very light exudate. If the system is not easily removed, then moisten or soak to assist removal.

Storage

Store at room temperature.

Antimicrobial Barrier Testing

Table 1: TheraBond[®]3D Antimicrobial Barrier System: Microbial Strike Through*

Challenge Organism	Inoculum Concentration Applied to Exterior Dressing Surface (Day 1)	Viable Organisms Recovered from Patient-Contact Dressing Surface (Day 14)
Staphylococcus aureus	2.44 x 10 ⁵	None (no growth)
Methicillin-resistant Staphylococcus aureus (MRSA)	2.19 x 10 ⁵	None (no growth)
Vancomycin-resistant E. faecalis (VRE)	1.51 x 10 ⁵	None (no growth)
Pseudomonas aeruginosa	1.01 x 10 ⁵	None (no growth)
Eschericia coli	1.04 x 10 ⁵	None (no growth)
Aspergillus brasiliensis	2.55 x 10 ⁵	None (no growth)
Candida albicans	2.40 x 10 ⁵	None (no growth)

* Data on file at Alliqua BioMedical, Inc.

Caution: Federal law restricts sale of the device to or on the order of a physician or licensed practitioner.

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* Datos archivados en Alliqua BioMedical, Inc.